

**General remarks**

The Cytometry Facility locations in room LLF D4.009 and LLF D4.020 are **BSL-2** working areas. No food and drinks inside the rooms!

All waste containers have to be decontaminated before disposing of any liquids.

**Solid Waste**

- Solid BSL-1 and BSL-2 waste is collected in the red waste container
- Do not fill the red waste containers with non-biohazardous waste
- Non-biohazardous waste can be disposed in the standard trash bins.  
Any trash placed in standard trash bins has to be safe for the cleaning personnel.

**Aria cell sorter (liquid waste)**

- **Decontaminated waste containers must be treated by the facility staff before disposal!**

***Last user of the day at Aria cell sorter***

- Disconnect full waste container from the instrument and attach an empty waste container to the instrument
- Add the appropriate volume of Bleach/Javel, (located under the sink) to the tank
- Put the magnetic “bleach added” label on the tank and incubate overnight.

**Cytek Aurora (liquid waste)**

- **Decontaminated waste containers must be treated by the facility staff before disposal!**

***User changing the waste tank when prompted by the software***

- Transfer the full waste container to the sink and add **200 ml** of Na-hypochlorite/Bleach/Javel (blue container located under the sink).
- Write a note on the tank that bleach has been added including date.

**Amounts of Na-hypochlorite/Bleach/Javel to add:**

Tank	Fill level	Amount of bleach to add
10 l plastic tanks (BD instruments)	¼ full	25 ml
	½ full	50 ml
	full	100 ml
20 l carton tanks (Aurora)	full	200 ml